



# Olympus™ Filter and Disinfectant Change Log

Freq.	Disinfectant Solution Test			Disinfectant Filter		3/1 micron "Pre-Sediment Filter"	.45 micron "Sediment Filter"	0.2 micron Internal Filter	Air Filter		Initials	Date
	A	B	C	Station A	Station B*				Station A	Station B		
	Daily			Every Disinfectant Change	Every 45 days (Depends on Water quality)	Every 45 days (Depends on Water quality)	Every 45 days (Depends on Water quality)	Every 3 months				
1	P F	P F	P F									
2	P F	P F	P F									
3	P F	P F	P F									
4	P F	P F	P F									
5	P F	P F	P F									
6	P F	P F	P F									
7	P F	P F	P F									
8	P F	P F	P F									
9	P F	P F	P F									
10	P F	P F	P F									
11	P F	P F	P F									
12	P F	P F	P F									
13	P F	P F	P F									
15	P F	P F	P F									
16	P F	P F	P F									
17	P F	P F	P F									
18	P F	P F	P F									
19	P F	P F	P F									
20	P F	P F	P F									
21	P F	P F	P F									
22	P F	P F	P F									
23	P F	P F	P F									
24	P F	P F	P F									
25	P F	P F	P F									
26	P F	P F	P F									
27	P F	P F	P F									
28	P F	P F	P F									
29	P F	P F	P F									
30	P F	P F	P F									
31	P F	P F	P F									

\* for dual basin washers